

ProCare Dent 11 A - 10 kg  
Powder cleaning detergent, alkaline, 10 kg  
for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for handpieces
- Adapted to the process sequences in Miele washer-disinfectors
- Excellent material compatibility
- Free from oxidative additives
- Corrosion protection for stainless steel instruments and trays

EAN: 4029577158217 / Material number: 12786540 / Old Material Number: 66911201V6

<b>Product category</b>	
Thermal disinfectors	•
<b>Applications</b>	
Removal of organic residues	•
<b>Dispensing</b>	
pH value at 20°C (68°F)	12,0
<b>Dispensing recommendation</b>	
Dispensing recommendation per rinse	3 g/l (30 g/cycle)
Dispensing recommendation (per cycle)	3 g/l (30 g/cycle)
Intensity of soiling	Strong
Suitable for all levels of water hardness	•
<b>Main applications</b>	
Dental practices	•
<b>Product characteristics</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Phosphate, Non-ionic surfactants, Silicates
Minimum storage temperature in °C	0
Maximum storage temperature in °C	30
<b>Product characteristics</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Phosphate, Non-ionic surfactants, Silicates
Minimum storage temperature in °C	0
Maximum storage temperature in °C	30
Suitable for standard instruments	•
Suitable for automatic reprocessing of instruments and utensils	•
High degree of material protection against corrosion	•
Good cleaning capability	•
Not suitable for handpieces	Suitable for standard instruments, Not suitable for transmission instr., For machine reprocessing of instruments
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	10 9/16 (268)
External dimensions, net width in in. (mm)	10 7/16 (265)
External dimensions, net depth in in. (mm)	10 7/16 (265)
Load size in g	352.73 (10000)
Gross weight in lb (kg)	0.35 (10)
<b>Compliance, test certificates, and symbols</b>	
CE	•
GHS / CLP compliant	•