

ProCare Dent 30 C - 5 l  
Neutralizing agent, acidic, 5 l  
for the automatic reprocessing of dental instruments and handpieces.



- Particularly suitable for transmission instruments
- Adapted to the process sequences in Miele washer-disinfectors
- Good material compatibility, free from phosphates
- Neutralizes alkali residues that are carried over
- Based on citric acid

EAN: 4029577158224 / Material number: 12796330 /  
Old Material Number: 66915303V6

Product category	
Thermal disinfectors	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
pH value at 20°C (68°F)	1,0
Dispensing recommendation	
Dispensing recommendation per rinse	1 ml/l (0.1%)
Dispensing recommendation (per cycle)	1 ml/l (0.1%)
Large opening for suction lance	•
Suitable for all levels of water hardness	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
Minimum storage temperature in °C	-3
Maximum storage temperature in °C	30
Dimensions and weight	
External dimensions, gross depth in mm	150.00
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
Minimum storage temperature in °C	-3
Maximum storage temperature in °C	30
Suitable for anodized aluminum	•
Suitable for standard instruments and handpieces	•
Does not contain surfactants	•
Very good material compatibility	•
Suitable for sensitive materials	•
Neutralization of alkaline carryover	•
Excellent rinsing properties	•
Does not contain phosphorus compounds	•
Prevention of limescale in wash cabinet	•
Not suitable for handpieces	Suitable for standard / transmission instruments, Suitable for anodized aluminum, Acid resistance Load items and drains
Dimensions and weight	
External dimensions, net height in in. (mm)	8 9/16 (217)
External dimensions, net width in in. (mm)	5 5/8 (142)
External dimensions, net depth in in. (mm)	9 13/16 (248)
External dimensions, gross height in in. (mm)	9 7/8 (250)
External dimensions, gross width in in. (mm)	7 1/2 (190)
External dimensions, gross depth in in. (mm)	5 15/16 (150)

ProCare Dent 30 C - 5 l  
Neutralizing agent, acidic, 5 l  
for the automatic reprocessing of dental instruments and handpieces.



EAN: 4029577158224 / Material number: 12796330 / Old Material Number: 66915303V6	
Gross weight in lb (kg)	0.21 (6)
Fill volume in fl.oz. (ml)	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•