

ProCare Dent 30 P - 5 l  
Neutralizing agent, acidic, 5 l  
for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for handpieces
- Adapted to the process sequences in Miele washer-disinfectors
- Removes discolouration and soiling which is soluble in acid
- Neutralizes alkali residues that are carried over
- Based on phosphoric acid

EAN: 4029577158231 / Material number: 12796380 / Old Material Number: 66914303V6

Product category	
Thermal disinfectors	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
pH value at 20°C (68°F)	0,7
Dispensing recommendation	
Dispensing recommendation per rinse	1 ml/l (0.1%)
Dispensing recommendation (per cycle)	1 ml/l (0.1%)
Large opening for suction lance	•
Suitable for all levels of water hardness	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphoric acid, Citric acid
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphoric acid,Citric acid
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
Suitable for porcelain / ceramics	•
Not suitable for anodized aluminum	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for light and non-ferrous metals	•
Suitable for standard instruments	•
Protects against corrosion	•
Does not contain surfactants	•
Neutralization of alkaline carryover	•
Excellent rinsing properties	•
Prevention of limescale in wash cabinet	•
Not suitable for handpieces	Suitable for standard instruments, Not suitable for light and non-ferrous metals, Not suitable for chrome-plated parts, Not suitable for nickel-plated parts, Not suitable for anodized aluminum, Not suitable for transmission instr., Suitable for porcelain/ceramics
Dimensions and weight	
External dimensions, net height in in. (mm)	8 9/16 (217)
External dimensions, net width in in. (mm)	5 5/8 (142)

ProCare Dent 30 P - 5 l  
Neutralizing agent, acidic, 5 l  
for the automatic reprocessing of dental instruments.



EAN: 4029577158231 / Material number: 12796380 / Old Material Number: 66914303V6

External dimensions, net depth in in. (mm)	9 13/16 (248)
Gross weight in lb (kg)	0.28 (8)
Fill volume in fl.oz. (ml)	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•