

## ProCare Dent 40 - 1 l [Typ 1]

Rinsing agent, 1 l

for the automatic reprocessing of dental instruments and handpieces.



- Suitable for dental instruments and handpieces
- Adapted to the process sequences in Miele washer-disinfectors
- Improves and accelerates drying
- Neutralizes alkali residues that are carried over
- Easier to add to the door dispenser thanks to the pouring aid

EAN: 4029577158255 / Material number: 12796420 / Old Material Number: 66917103V6

Product category	
Thermal disinfectors	•
Dispensing	
pH value at 20°C (68°F)	2,0
Dispensing recommendation	
Dispensing recommendation per rinse	0.3 ml/l (0.03%)
Dispensing recommendation (per cycle)	0.3 ml/l (0.03%)
Suitable for all levels of water hardness	•
Small opening as pouring aid	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Rinsing agent
Aggregate form	Liquid concentrate
Ingredients	Preservatives, Surfactants, Citric acid
Minimum storage temperature in °C	0
Maximum storage temperature in °C	30
Dimensions and weight	
External dimensions, gross depth in mm	100.00
Product characteristics	
Medium type	Rinsing agent
Aggregate form	Liquid concentrate
Ingredients	Preservatives, Surfactants, Citric acid
Minimum storage temperature in °C	0
Maximum storage temperature in °C	30
Suitable for porcelain / ceramics	•
Suitable for standard instruments and handpieces	•
Suitable for all metals used in medical devices	•
Improves and shortens the drying process	•
Neutralization of alkaline carryover	•
Not suitable for handpieces	Suitable for standard / transmission instruments, Suitable for porcelain/ceramics, For all metals in medical devices
Prevention of water stains	•
Dimensions and weight	
External dimensions, net height in in. (mm)	9 1/2 (241)
External dimensions, net width in in. (mm)	4 3/16 (106)
External dimensions, net depth in in. (mm)	2 5/16 (58)
External dimensions, gross height in in. (mm)	10 1/4 (260)
External dimensions, gross width in in. (mm)	2 3/8 (60)
External dimensions, gross depth in in. (mm)	3 15/16 (100)
Gross weight in lb (kg)	0.04 (1)
Fill volume in fl.oz. (ml)	33.81 (1000) fl oz (US) (Milliliter)

ProCare Dent 40 - 1 l [Typ 1]  
Rinsing agent, 1 l  
for the automatic reprocessing of dental instruments and handpieces.



EAN: 4029577158255 / Material number: 12796420 / Old Material Number: 66917103V6

Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•