

E 315/1

Full demineralization single use resin

For the regeneration of full demineralization cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



- For use in the LP 2800/1 empty cartridge
- 20 L of homogeneously mixed resins
- Carton with 2 x 10 L bags
- Vacuum sealed in plastic bags
- Including replacement filter bag

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•

E 315/1

Full demineralization single use resin

For the regeneration of full demineralization cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PLW 8683	•
PLW 8683 CD	•
PWD 8682	•
PWD 8682 C	•
PWD 8692	•
Product category	
Water treatment	•
Dimensions and weight	
External dimensions, gross height in in. (mm)	11 1/4 (285)
External dimensions, gross width in in. (mm)	10 1/4 (260)
External dimensions, gross depth in in. (mm)	17 3/8 (440)
Net weight in lb (kg)	30 (13.77)
Gross weight in lb (kg)	31 (14.27)