

ProCare Lab 30 P - 5 l  
Neutralizing agent based on phosphoric acid, acidic, 5 L  
for machine cleaning of laboratory glassware and utensils.



- To be used as neutralizing agent or acid pre-cleaning product
- Does not contain surfactants
- Specially developed for Miele Professional laboratory glasswashers
- Suitable for stainless steel, laboratory glassware, ceramics, and plastics
- Removes stubborn residues that are soluble in acids

EAN: 4029577160364 / Material number: 12960740 /  
Old Material Number: 62928102V2

<b>Product category</b>	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
<b>Applications</b>	
Removal of limescale and other acid-soluble impurities	•
<b>Dispensing</b>	
Recommended dosage per wash cycle	1-4 ml/l (0.1-0.4%)
Large opening for suction lance	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
<b>Main applications</b>	
Hospital laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
<b>Product characteristics</b>	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphate, Phosphoric acid, Citric acid
pH value at 20°C (68°F)	0,7
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for aluminum	•
Not suitable for light and non-ferrous metals	•
Not suitable for light alloys	•
Does not contain surfactants	•
Free from microplastics	•
Acidic pre-cleaning	•
Excellent rinsing properties	•
Neutralization of alkaline carryover	•
<b>Storage information</b>	
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
<b>Dimensions and weight</b>	
External dimensions, net height in mm	8 9/16 (217)
External dimensions, net width in mm	5 5/8 (142)
External dimensions, net depth in mm	9 13/16 (248)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
<b>Compliance, test certificates, and symbols</b>	
GHS / CLP compliant	•

ProCare Lab 30 P - 5 l

Neutralizing agent based on phosphoric acid, acidic, 5 L  
for machine cleaning of laboratory glassware and utensils.