

E 338

Insert

for the optimum loading of up to 8 half-shell trays.

EAN: 4002511116160 / Material number: 03080600 / Old Material Number: 69533801D



- Suitable for 8 half-shell trays
- equipped with 10 holders (8 compartments), W 295, D 33 mm
- Maximum tray size 290 x 30 mm
- For use in upper and lower basket

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	8
<b>Capacity</b>	
Trays [quantity]	8
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	110
External dimensions, net width in mm	305
External dimensions, net depth in mm	445
External dimensions, gross height in mm	125
External dimensions, gross width in mm	310
External dimensions, gross depth in mm	500
Net weight in kg	0,62
Gross weight in kg	0,84