


O 183

Upper basket

for the optimum loading of instruments for MIS, arthroscopy and urology.



- Left side: 10 injector nozzles or free accommodation of inserts
- Right side: 14 injector nozzles/ funnels for hollow instruments
- Vertical clearance 285 +/- 30mm; hollow instruments 370 +/- 30mm
- Powder dispensing not possible

EAN: 4002511114920 / Material number: 03807470 / Old Material Number: 69118301D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Capacity	
Lumened instruments [number]	14
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	270
External dimensions, net width in mm	475
External dimensions, net depth in mm	530
External dimensions, gross height in mm	410
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	555
Net weight in kg	4,00
Gross weight in kg	4,80