


E 126

Insert

For the optimum loading of infusion flasks. Surface area 45 x 45 mm.



- 48 compartments, dimension: 45 x 45 mm (bottle), 28 x 28 mm (neck)

EAN: 4002511115835 / Material number: 03808520 / Old Material Number: 69512601D

Product classification	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	48
Area of application	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	83
External dimensions, net width in mm	445
External dimensions, net depth in mm	445
External dimensions, gross height in mm	160
External dimensions, gross width in mm	475
External dimensions, gross depth in mm	475
Net weight in kg	1,92
Gross weight in kg	2,05