

E 129

Insert

For the optimum loading of infusion flasks. Total area 84 x 84 mm.



- 20 compartments, dimension: 84 x 84 mm (bottle), 46 x 46 mm (neck)

EAN: 4002511115866 / Material number: 03808550 / Old Material Number: 69512901D

| Product classification | |
|---|-----------------|
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| Product category | |
| Insert for various flasks | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Number of compartments | 20 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 113 |
| External dimensions, net width in mm | 445 |
| External dimensions, net depth in mm | 445 |
| External dimensions, gross height in mm | 480 |
| External dimensions, gross width in mm | 480 |
| External dimensions, gross depth in mm | 130 |
| Net weight in kg | 1,41 |
| Gross weight in kg | 1,78 |