

E 134  
Insert  
for optimum loading of 210 slides.



- suitable for 210 slides
- Equipped with 26 compartments at 11 mm x 26 mm, wire thickness 3 mm

EAN: 4002511115897 / Material number: 03808600 / Old Material Number: 69513401D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8582 CD	•
PG 8504	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	26
Compartment width in mm	11
Compartment length in mm	26
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Slides [number]	210
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	73
External dimensions, net width in mm	200
External dimensions, net depth in mm	445
External dimensions, gross height in mm	80
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	460
Net weight in kg	1,50
Gross weight in kg	1,61