

E 135/2

Insert

for optimum loading of 19 baby bottles up to 90 ml.

EAN: 4002511115927 / Material number: 03808650 / Old Material Number: 69513503D



- Suitable for 19 baby bottles x 90 ml
- Suitable for bottles sized 51 x 51 mm, neck: 42 x 42 mm
- For upper basket A 102 and lower basket A 151

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for various flasks	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of baby bottles	•
Capacity	
Baby bottles [number]	19
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	125
External dimensions, net width in mm	192
External dimensions, net depth in mm	447
External dimensions, gross height in mm	140
External dimensions, gross width in mm	235
External dimensions, gross depth in mm	480
Net weight in kg	1,97
Gross weight in kg	2,10