

E 331
Mobile unit
for the optimum loading of butyrometers.



- equipped with 39 injector nozzles

EAN: 4002511116115 / Material number: 03809350 / Old Material Number: 69533101D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Other mobile units	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of butyrometers	•
Capacity	
Butyrometers [number]	39
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	517
External dimensions, gross height in mm	500
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	530
Net weight in kg	5,20
Gross weight in kg	7,10
Companion pack	
SD-B injector nozzle for butyrometer [number]	39
Rinse nozzle for powder dispenser [number]	1