

E 351

Injector nozzle

For the optimum cleaning of the wash items. Size 4 x 160 mm.



- Suitable for use in a mobile injector unit
- Can be combined with clip E 353
- Ø 4 x 160 mm, screw connection

EAN: 4002511117075 / Material number: 03809500 / Old Material Number: 69735101D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	4
External dimensions, net depth in mm	4
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	4
Net weight in kg	0,01
Gross weight in kg	0,02