

E 354

Clamps

Height-adjustable for optimum positioning of laboratory glassware on the nozzle.

EAN: 4002511117105 / Material number: 03809540 / Old Material Number: 69735401D



- Suitable for injector nozzle E 352
- Height adjustable locking
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

| | |
|--|---|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| Product category | |
| Other accessories | • |
| Product properties | |
| Material | Wire, rilsan-coated, Stainless steel |
| Colour | Stainless Steel, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0,00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 150 |
| External dimensions, net width in mm | 70 |
| External dimensions, net depth in mm | 70 |
| External dimensions, gross height in mm | 160 |
| External dimensions, gross width in mm | 100 |
| External dimensions, gross depth in mm | 150 |
| Net weight in kg | 0,05 |
| Gross weight in kg | 0,05 |