

E 364

Insert

for the optimum loading of 36 wide-necked teats.



- Suitable for 36 wide-necked teats
- Equipped with 36 compartments at 44 mm x 44 mm
- Flip-up and lockable lid

EAN: 4002511116290 / Material number: 03809640 / Old Material Number: 69536401D

<b>Product classification</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
Compartment width in mm	44
Compartment length in mm	44
<b>Area of application</b>	
Reprocessing of baby bottles	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	77
External dimensions, net width in mm	215
External dimensions, net depth in mm	445
External dimensions, gross height in mm	80
External dimensions, gross width in mm	235
External dimensions, gross depth in mm	460
Net weight in kg	1,48
Gross weight in kg	1,91