

E 381

Mobile unit

for the optimum loading of intubation material.



- 20 nozzles Ø 4.0 x 30 mm, 5 nozzles Ø 2.5 x 30 mm (5 loose)
- 6 x E 466 injector nozzles for breathing bags Ø 8.0 x 333 mm
- 2 x E 431 injector nozzles for bellows, Ø 8.0 x 193 mm
- Equipped with connection for hot air drying
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002511116474 / Material number: 03809890 / Old Material Number: 69538101D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Mobile unit for anaesthetic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	502
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	540
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	500
Net weight in kg	6,80
Gross weight in kg	8,70