

E 137

Insert

for optimum loading of 56 half petri dishes.



- Suitable for 56 half petri dishes with Ø 100 mm
- Equipped with 56 holders, at 70 mm height, spacing approx. 26 mm
- For use in lower basket
- For use with insert E 136

EAN: 4002511115958 / Material number: 03830290 / Old Material Number: 69513701D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for half petri dishes	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	26
Rail height in mm	70
Capacity	
Half petri dishes [number]	56
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0

E 137

Insert

for optimum loading of 56 half petri dishes.

EAN: 4002511115958 / Material number: 03830290 / Old Material Number: 69513701D



Dimensions and weight	
External dimensions, net height in mm	485
External dimensions, net width in mm	445
External dimensions, net depth in mm	95
External dimensions, gross height in mm	140
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	580
Net weight in kg	2,03
Gross weight in kg	3,45