

E 443

Irrigation sleeves

for the optimum cleaning of MIS instruments.

EAN: 4002511908192 / Material number: 04203840 / Old Material Number: 69744301D



- Suitable for MIS instruments with \varnothing 8 – 8.5 mm, screw connection
- Length 121 mm mm, \varnothing 10 mm mm

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Irrigation sleeve	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	121
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Diameter in mm	10
Net weight in kg	0,03
Gross weight in kg	0,03