

E 800  
Insert  
for optimum loading of 3 mesh or kidney trays.



- Equipped with 4 holders
- For use in upper and lower basket

EAN: 4002512087773 / Material number: 04637720 / Old Material Number: 69580001D

<b>Product classification</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for mesh trays	•
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of kidney dishes	•
<b>Capacity</b>	
Kidney dishes [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	165
External dimensions, net width in mm	140
External dimensions, net depth in mm	290
External dimensions, gross height in mm	190
External dimensions, gross width in mm	195
External dimensions, gross depth in mm	430
Net weight in kg	0,29
Gross weight in kg	0,49