

DOS S

Dispenser pump rack

With slanting base for the complete emptying of the liquid media containers.

EAN: 4002512161107 / Material number: 04814240 / Old Material Number: 59170031D



- Suitable for the installation of the dispenser pump casing DOS G2
- Suitable for the installation of the dispenser pump unit DOS GP3
- With four swivel castors for maximum mobility

Product classification	
Washing machines	•
Machine/accessory category	
PW 413	•
PW 418	•
PW 5064	•
PW 5062	•
PW 5055	•
PW 5065	•
PW 5064 MOPSTAR 60	•
PW 5070	•
PW 5082	•
PW 5084 MOPSTAR 80	•
PW 5104 MOPSTAR 100	•
PW 5105 Vario	•
PW 5134 MOPSTAR 130	•
PW 6055	•
PW 6065	•
PW 6055 Vario	•
PW 6065 Vario	•
PW 6065 Plus	•
PW 6068 Plus	•
PW 6080	•
PW 6088	•
PW 6107	•
PW 6137	•
PW 5136	•
PW 6163	•
PW 6167	•
PW 6207	•
PW 6241	•
PW 6243	•
PW 6321	•
PW 6323	•
PW 811	•
PW 814	•
PW 818	•
PW 9055	•
PW 9065	•
PWT 6089	•
Product category	
Dispenser system	•
Product properties	
Electrical connection	OHNE
Material	Powder coated galvanised steel
Colour	Iron grey
Dimensions and weight	
External dimensions, net height in mm	1030
External dimensions, net width in mm	1000

DOS S

Dispenser pump rack

With slanting base for the complete emptying of the liquid media containers.

EAN: 4002512161107 / Material number: 04814240 / Old Material Number: 59170031D



External dimensions, net depth in mm	371
External dimensions, gross height in mm	1030
External dimensions, gross width in mm	1300
External dimensions, gross depth in mm	413
Net weight [kg]	9,00
Gross weight in kg	10,00