

E 466

Injector nozzle

for the optimum loading of breathing bags 8 x 333 mm.



EAN: 4002512223744 / Material number: 04966620 / Old Material Number: 69746601D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	333
External dimensions, gross height in mm	140
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	150
Diameter in mm	8
Net weight in kg	0,85
Gross weight in kg	1,00