

E 468

Insert

for the optimum loading of different utensils.



- Made of welded stainless steel mesh with lid and sub-divider
- Mesh width 5 x 5 x 1 mm

EAN: 4002512277280 / Material number: 05043640 / Old Material Number: 69546801D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	76
External dimensions, net width in mm	170
External dimensions, net depth in mm	250
External dimensions, gross height in mm	100
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	315
Net weight in kg	0,67
Gross weight in kg	0,82