

E 485  
Bracket  
for the optimum loading of 3 kidney trays.



EAN: 4002512398725 / Material number: 05109210 / Old Material Number: 69548501D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PG 8593	•
Product category	
Insert for kidney dishes/basins	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of kidney dishes	•
Capacity	
Kidney dishes [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	195
External dimensions, net width in mm	80
External dimensions, net depth in mm	464
External dimensions, gross height in mm	70
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	505
Net weight in kg	0,41
Gross weight in kg	0,56