

E 486  
Bracket  
for the optimum loading of basins.



- Suitable for 1 basins

EAN: 4002512398732 / Material number: 05109230 / Old Material Number: 69548601D

<b>Product classification</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for kidney dishes/basins	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of basins	•
<b>Capacity</b>	
Basins [number]	1
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	275
External dimensions, net width in mm	464
External dimensions, net depth in mm	10
External dimensions, gross height in mm	475
External dimensions, gross width in mm	465
External dimensions, gross depth in mm	25
Net weight in kg	0,21
Gross weight in kg	0,43