


AUF 1 Holder

for optimum reprocessing of transmission instruments in upper basket O 177/1.



- Comprises: seal, insert nut, clamp
- Rack lower part, rack upper part, 5 filter plates

EAN: 4002512427517 / Material number: 05159510 / Old Material Number: 68551101D

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8591	•
Product category	
Holder and/or adapter for transmission instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	80
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	100
Net weight in kg	0,90
Gross weight in kg	1,84
Companion pack	
FP filter plate	•