

E 730  
 Insert  
 for the optimum loading of 20 theatre shoes.



EAN: 4002512423120 / Material number: 05271020 / Old Material Number: 69573001D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	20
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	290
External dimensions, net width in mm	470
External dimensions, net depth in mm	540
External dimensions, gross height in mm	420
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	600
Net weight in kg	8,70
Gross weight in kg	10,00