


E 489

Universal bracket

for loading diverse wash items such as insoles for theatre shoes.



EAN: 4002512431835 / Material number: 05337120 / Old Material Number: 69548901D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	60
External dimensions, net width in mm	464
External dimensions, net depth in mm	10
External dimensions, gross height in mm	35
External dimensions, gross width in mm	130
External dimensions, gross depth in mm	470
Net weight in kg	0,12
Gross weight in kg	0,23