

E 985  
Bracket  
for optimum loading of 2 kidney trays.



- For use in insert E 984

EAN: 4002512431330 / Material number: 05337760 / Old Material Number: 69598501D

Product classification	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
Product category	
Insert for kidney dishes/basins	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of kidney dishes	•
Capacity	
Kidney dishes [number]	2
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	70
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	505
Net weight in kg	0,23
Gross weight in kg	0,43