

E 986  
Bracket  
for the optimum loading of 5 basins.



- For use in insert E 984

EAN: 4002512431347 / Material number: 05337770 / Old Material Number: 69598601D

| Product classification                      |                 |
|---|-----------------|
| Large-capacity washer-disinfectors, medical | •               |
| Machine affiliation                         |                 |
| PG 8527                                     | •               |
| PG 8528                                     | •               |
| Product category                            |                 |
| Insert for kidney dishes/basins             | •               |
| Product properties                          |                 |
| Material                                    | Stainless steel |
| Colour                                      | Stainless Steel |
| Area of application                         |                 |
| Reprocessing of basins                      | •               |
| Capacity                                    |                 |
| Basins [number]                             | 1               |
| Standard electrical connection              |                 |
| Number of phases                            | 0               |
| Voltage in V                                | 0               |
| Frequency in Hz                             | 0               |
| Total rated load in kW                      | 0,00            |
| Fuse rating in A                            | 0               |
| Dimensions and weight                       |                 |
| External dimensions, net height in mm       | 0               |
| External dimensions, net width in mm        | 0               |
| External dimensions, net depth in mm        | 0               |
| External dimensions, gross height in mm     | 25              |
| External dimensions, gross width in mm      | 330             |
| External dimensions, gross depth in mm      | 435             |
| Net weight in kg                            | 0,19            |
| Gross weight in kg                          | 0,36            |