

E 987

Bracket

for optimum loading of 15 theatre shoes up to size 44.



- For use in insert E 984

EAN: 4002512431354 / Material number: 05337800 / Old Material Number: 69598701D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	15
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	30
External dimensions, gross width in mm	460
External dimensions, gross depth in mm	475
Net weight in kg	0,23
Gross weight in kg	0,40