


E 989

Universal bracket

for 8 insoles for theatre shoes, for example.



- For use in insert E 984

EAN: 4002512431828 / Material number: 05376740 / Old Material Number: 69598901D

<b>Product classification</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8528	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Insoles [number]	8
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	20
External dimensions, gross width in mm	330
External dimensions, gross depth in mm	415
Net weight in kg	0,10
Gross weight in kg	0,26