


A 5

Lid

For optimum securing of delicate items.



- Stainl. steel frames w/ plastic netting (material quality PA 12)
- For use in O 184 and U 184

EAN: 4002513126310 / Material number: 05637190 / Old Material Number: 69500501D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8535	•
PG 8536	•
Product category	
Cover	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated, Plastic
Colour	White
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	8
External dimensions, net width in mm	280
External dimensions, net depth in mm	280
External dimensions, gross height in mm	20
External dimensions, gross width in mm	470
External dimensions, gross depth in mm	470
Net weight in kg	0,40
Gross weight in kg	0,55