

ML/2

Magnetic strip

for automatic recognition of mobile unit/spray arm monitoring



- ML/2 is required for spray arm monitoring
- 5 magnets, individually configurable
- 15 possible combinations

EAN: 4002513274158 / Material number: 05723270 / Old Material Number: 69745703D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Magnetic strip for automatic mobile unit recognition	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel, Black
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	9
External dimensions, net width in mm	25
External dimensions, net depth in mm	275
External dimensions, gross height in mm	95
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	315
Net weight in kg	0,17
Gross weight in kg	0,32