

## DOS G 10

Dispenser module for 5-10 litre canisters

With level indicator for the automatic dispensing of liquid cleaning agents.



- For dispensing rinsing agent
- 100 mm siphon for 5-10 l containers
- Conversion kit for siphons for 10-30 l canisters available
- Dispenses up to 100 ml/60 seconds

EAN: 4002513283259 / Material number: 05729360 / Old Material Number: 69747006D

<b>Product classification</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
<b>Product category</b>	
Dispenser system	•
<b>Standard electrical connection</b>	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	220-240
Frequency in Hz	50-60
Total rated load in kW	0,02
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	118
External dimensions, net width in mm	235
External dimensions, net depth in mm	115
External dimensions, gross height in mm	310
External dimensions, gross width in mm	420
External dimensions, gross depth in mm	250
Net weight in kg	2,42
Gross weight in kg	2,81