

E 492  
Insert  
for the optimum loading of kidney trays.



- equipped with 9 holders, height 86 mm mm, distance 49 mm mm
- For use in lower basket

EAN: 4002513317343 / Material number: 05771140 / Old Material Number: 69549201D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for kidney dishes/basins	•
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	49
Rail height in mm	86
<b>Area of application</b>	
Reprocessing of kidney dishes	•
<b>Capacity</b>	
Kidney dishes [number]	9
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, net width in mm	256
External dimensions, net depth in mm	474
External dimensions, gross height in mm	125
External dimensions, gross width in mm	320
External dimensions, gross depth in mm	480
Net weight in kg	1,06
Gross weight in kg	1,31