


UC 30-90/70-78

Box plinth

to raise the machine, 70 cm deep, 90 cm wide.



- Suitable for all machines, width 90, depth 70 cm
- Can be screwed onto the appliance

EAN: 4002513347807 / Material number: 05818100 / Old Material Number: 69379009D

Product classification	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
Product category	
Plinth/base	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	300
External dimensions, net width in mm	900
External dimensions, net depth in mm	700
External dimensions, gross height in mm	370
External dimensions, gross width in mm	1000
External dimensions, gross depth in mm	730
Net weight in kg	13,20
Gross weight in kg	14,40