

U 168/1  
Lower basket  
for the ideal loading of theatre shoes.



- Suitable for 20 theatre shoes up to size 45
- 20 holders 295 mm mm
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002513390940 / Material number: 05878590 / Old Material Number: 69116803D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Insert for theatre shoes	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Rail height in mm	295
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	20
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	315
External dimensions, net width in mm	535
External dimensions, net depth in mm	515
External dimensions, gross height in mm	320
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	545
Net weight in kg	3,20
Gross weight in kg	5,20