


E 902/1

Mobile unit

for loading up to 3 MIS theatre sets in a modular system.



- Suitable for 6 insert modules with integrated nozzles/adapters
- Upper & lower levels free for inserts, e.g. E 142
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets ML)
- Simultaneous cleaning of mesh trays and MIS theatre sets

EAN: 4002513505542 / Material number: 06059530 / Old Material Number: 69590202D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for MIS instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Capacity	
DIN mesh trays [number]	3
SPRI mesh trays [number]	1
MIS sets [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	620
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	650
External dimensions, gross width in mm	650
External dimensions, gross depth in mm	800
Net weight in kg	24,00
Gross weight in kg	34,00
Companion pack	
E 362 blind stopper [number]	6
E 444 drum for fibre optic cables/suction hoses [number]	3
E 445 – 12 covers for E 442 irrigation sleeve, opening with Ø 6 mm [number]	2
E 446 – 12 covers for E 443 irrigation sleeve, opening with Ø 10 mm [number]	1
E 905/1 module insert for short MIS instruments [number]	3
E 906/1 module insert for long MIS instruments [number]	3
2 m silicone hose, Ø 5 mm	•
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1