


U 175/1

Lower basket

for the optimum loading of narrow necked glassware.



- Equipped with 33 injector nozzles
- Vertical clearance: 170 mm mm\*
- Only suitable for use with upper baskets O 175 / O 187 / O 184

\* Depending on the laboratory glassware the max. vertical clearance can be lower.

EAN: 4002513805543 / Material number: 06439170 / Old Material Number: 69117503D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	170
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	33
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	475
External dimensions, net depth in mm	531
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	6,40
Gross weight in kg	7,65
<b>Companion pack</b>	
E 351 injector nozzle, Ø 4 mm, length 160 mm [number]	33
E 353 clip for injector nozzle E 351 [number]	33
Rinse nozzle for powder dispenser [number]	1