

E 715/1

Mobile unit

for holding up to 4 AN sets and/or 12 breathing tubes (1.5 m) and accessories.



- Holds breathing tubes on a spiral rack
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513872699 / Material number: 06491890 / Old Material Number: 69571502D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Mobile unit for anaesthetic instruments	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Capacity</b>	
AN sets [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	530
External dimensions, net width in mm	600
External dimensions, net depth in mm	645
External dimensions, gross height in mm	560
External dimensions, gross width in mm	640
External dimensions, gross depth in mm	680
Net weight in kg	9,70
Gross weight in kg	10,00
<b>Companion pack</b>	
E 430/1 1/3 mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	3
E 434 holder for paediatric breathing tubes [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	6
E 496 injector nozzle 4 x 120 mm [number]	15
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1