

E 479

Holder

for micro instruments up to 4 mm Ø.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513918175 / Material number: 06553750 / Old Material Number: 69747901D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Holder	•
<b>Product properties</b>	
Material	Plastic
Colour	Transparent
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	5
External dimensions, net width in mm	21
External dimensions, net depth in mm	21
External dimensions, gross height in mm	95
External dimensions, gross width in mm	110
External dimensions, gross depth in mm	80
Net weight in kg	0,19
Gross weight in kg	0,23