

E 494  
 Insert  
 for the gentle holding of 5 microtiter plates.



EAN: 4002513926064 / Material number: 06570920 / Old Material Number: 69549401D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
<b>Capacity</b>	
Microtiter plates [number]	5
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	35
External dimensions, net width in mm	205
External dimensions, net depth in mm	440
External dimensions, gross height in mm	100
External dimensions, gross width in mm	265
External dimensions, gross depth in mm	180
Net weight in kg	0,85
Gross weight in kg	1,00