


E 522/1

Insert

for the optimum loading of 9 impression trays.



- Equipped with 9 holders

EAN: 4002514072722 / Material number: 06697010 / Old Material Number: 68552202D

<b>Product classification</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Impression trays [number]	9
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	140
External dimensions, net width in mm	100
External dimensions, net depth in mm	190
External dimensions, gross height in mm	170
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	240
Net weight in kg	0,26
Gross weight in kg	0,36