

DOS K 60

Dispenser module for 5-10 litre canisters

With level indicator for the automatic dispensing of liquid cleaning agents.

EAN: 4002514254890 / Material number: 06940330 / Old Material Number: 69747008D



- Dispenser monitoring – in accordance with DIN EN ISO 15883
- For dispensing liquid cleaning agent/chemical disinfectant
- 300 mm siphon for 5-10 litre containers
- Conversion kit for siphons for 10-30 l canisters available
- Dispenses up to 120 Milliliter ml/60 seconds

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
Product category	
Dispenser system	•
Product properties	
Production rate in ml/60 s	120
Length of suction lance in mm	333
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	230-240
Frequency in Hz	50
Total rated load in kW	0,02
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	118
External dimensions, net width in mm	235
External dimensions, net depth in mm	115
External dimensions, gross height in mm	290
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	420
Net weight in kg	2,45
Gross weight in kg	2,70