

DOS K 60/1

Dispenser module for 5 litre canister

With level indicator for the automatic dispensing of liquid cleaning agents.

EAN: 4002514254906 / Material number: 06940340 / Old Material Number: 69747009D



- Dispenser monitoring – in accordance with DIN EN ISO 15883
- For dispensing liquid cleaning agent/chemical disinfectant
- 200 mm siphon for 5 litre containers
- Conversion kit for siphons for 10-30 l canisters available
- Dispenses up to 120 Milliliter ml/60 seconds

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
<b>Product category</b>	
Dispenser system	•
<b>Product properties</b>	
Production rate in ml/60 s	120
Length of suction lance in mm	200
<b>Standard electrical connection</b>	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	230-240
Frequency in Hz	50
Total rated load in kW	0,02
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	118
External dimensions, net width in mm	235
External dimensions, net depth in mm	115
External dimensions, gross height in mm	290
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	420
Net weight in kg	2,60
Gross weight in kg	3,00