

E 906/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.



- Holds up to 10 hollow instruments
- Subdivision for arthroscopy and laparoscopy, for example.

EAN: 4002514256313 / Material number: 06943250 / Old Material Number: 69590602D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
<b>Product category</b>	
Insert/module for MIS instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Capacity</b>	
Lumened instruments [number]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	40
External dimensions, net width in mm	461
External dimensions, net depth in mm	510
External dimensions, gross height in mm	160
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	530
Net weight in kg	2,87
Gross weight in kg	3,56
<b>Companion pack</b>	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	1
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	5
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	2
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
E 456 spring stay for MIS instruments [number]	2
E 464 holder for injector nozzle E 454 [number]	1
E 908/1 insert for separable MIS instruments/work inserts [number]	1