

E 955

Insert

for the optimum loading of drill sleeves with \varnothing 3 - 15 mm, length max. 125 mm



EAN: 4002514312057 / Material number: 07025330 / Old Material Number: 69595501D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of orthopaedic surgery instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	85
External dimensions, net width in mm	221
External dimensions, net depth in mm	156
External dimensions, gross height in mm	170
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	250
Net weight in kg	1,00
Gross weight in kg	1,00