

E 980  
 Adapter  
 for the optimum loading of long hollow instruments.



- For loading long hollow instruments
- Possible to screw in various nozzles or Luer lock adapters

EAN: 4002514312071 / Material number: 07025360 / Old Material Number: 69598001D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of lumened instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	15
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	80
Net weight in kg	0,04
Gross weight in kg	0,05