

## HUB 5001-M

Induction compatible gourmet oven dish

For frying, braising and gratinating.

EAN: 4002514463285 / Material number: 07247750 / Old Material Number: 22996230D



- Max. cap. approx 5 kg – Ext. size:  
WxDxH 44.8 x 26.8 x 8.6 cm
- Made from die-cast aluminium,  
hand-cast
- Easy to clean thanks to high-quality  
non-stick finish
- For ovens and all hobs except gas  
hobs
- Without stainless steel lid HBD  
60-22